

gRNA Synthesis Protocol

STEP 1: Find all 23bp genomic sites of the form 5'-N₂₀NGG-3' near your intended target site (ideally ±50bp). These may reside on the + or - strand.



STEP 2: Using NCBI blast, select sequences for which none or very few sequences of the form 5'-NNNNN NNNNN NNNNN NNNNN NGG-3' exist at any other location in the human genome (here the B's represent the actual bases at the target genomic location).

Option A

Option B

STEP 3: Incorporate 19bp of the selected target sequence as highlighted here: 5'-NNNNN NNNNN NNNNN NNNNN NGG-3' into the DNA fragment as indicated below:

```
TGTACAAAAAAGCAGGCTTTAAAGGAACCAATTCAGTCGACTGGATCCGGTACC
AAGGTCGGGCAGGAAGAGGGCCTATTTCCCATGATTCCCTTCATATTTGCATATA
CGATACAAGGCTGTTAGAGAGATAATTAGAATTAATTTGACTGTAAACACAAAG
ATATTAGTACAAAATACGTGACGTAGAAAAGTAATAATTTCTTGGGTAGTTTGA
GTTTTAAATATATTTTTAAATGGACTATCATATGCTTACCCTAACTTAAAG
TATTTCGATTTCTTGGCTTATATATCTTGTGGAAAGGACGAAACACCGNNNNN
NNNNNNNNNNNNNNGTTTTAGAGCTAGAAATAGCAAGTTAAATAAAGGCTAGTC
CGTTATCAACTTGAAAAGTGGCACCAGTCCGGTCTTTTTTCTAGACCCAGC
TTTTTCTAGACCCAGC
TTTTCTGTACAAAGTTGGCATTAA
```

Step 4: This 455bp fragment bears all components necessary for gRNA expression, namely: U6 promoter + target sequence + guide RNA scaffold + termination signal. Synthesize this as a gBlock from IDT (<http://www.idtdna.com/pages/products/genes/gblocks-gene-fragments>)

Step 5: Clone the synthesized gBlock into an empty backbone vector such as pCR-Blunt II-TOPO from Invitrogen (<http://products.invitrogen.com/ivgn/product/K280020>), or directly pcr amplify this fragment (gRNA_F: TGTACAAAAAAGCAGGCTTTAAAG, gRNA_R: TAATGCCAACCTTTGTACAAGAAAG) for transfection and gRNA expression.

STEP 3: Incorporate 19bp of the selected target sequence as highlighted here: 5'-NNNNN NNNNN NNNNN NNNNN NGG-3' into two 60mer oligonucleotides as indicated below (sequences are 5' to 3', and the regions marked in green and red are reverse complements of each other):

Insert_F:
TTTCTTGGCTTTATATATCTTGTGGAAAGGACGAAACACCGNNNNNNNNNNNN
NNNNNNNN
Insert_R:
GACTAGCCTATTTTAACTTGCTATTTCTAGCTCTAAACNNNNNNNNNNNN
NNNNNNNC

Step 4: Anneal the two oligos and extend these to make a 100bp double stranded DNA fragment using Phusion polymerase (<http://www.neb.com/nebecomm/products/productm0530.asp>).

Step 5: Linearize the gRNA cloning vector (<http://www.addgene.org/41824/>) using *AflI* and incorporate the 100bp DNA fragment from Step 4 above into it using Gibson assembly (<http://www.neb.com/nebecomm/products/productE2611.asp>). The resulting vector is the desired gRNA expression vector. **Note:** This synthesis strategy is amenable to construction of multiple gRNA expression vectors, and even large gRNA libraries using oligonucleotide pools synthesized using custom DNA arrays (refer Supplementary Fig. 11: <http://www.sciencemag.org/content/suppl/2013/01/03/science.1232033.DC1/Mali.SM.pdf>).

